SEE THE ENTIRE ERGODYNE[®] PROFLEX[®] SUPPORTS LINE OF BACK SUPPORTS // WRIST SUPPORTS // WRAPS & SLEEVES ON ERGODYNE.COM

ergodyne[®]

INTERGALACTIC HEADQUARTERS
1021 BANDANA BOULEVARD EAST // SUITE 220 // SAINT PAUL, MN 55108
TOLL-FREE 800 225 8238 // INT'L +1 651 642 9889 // ERGODYNE.COM

©2023 ERGODYNE. ALL RIGHTS RESERVED. ALL WRONGS REVERSED.
DESIGNED BY ERGODYNE IN MINNESOTA

CE / UKCA AUTHORIZED REPRESENTATIVE - REPRESENTED BY ARC:
AUTHORIZED REP COMPLIANCE LTD., GROUND FLOOR,
71 LOWER BAGGOT STREET, DUBLIN, D02 P593, IRELAND
AUTHORIZED REP COMPLIANCE LTD.,
ARC HOUSE, THURNHAM, LANCASTER, LA2 0DT, UK

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PROFLEX[®] BACK SUPPORT BRACE WITH COOLING/WARMING PACK

PROFLEX[®] BACK SUPPORT



HOW TO PUT ON YOUR PROFLEX[®] BACK SUPPORT:

1. Loosen suspenders and slip on the back support.
2. Grasp each end of the waistband. Stretch the right end with hook-and-loop fastener across the left and attach (approximately 4-6 inches (10-15 cm) of overlap). Model 1400UN universal back support may overlap more than 4-6 inches (10-15 cm) when worn on people with small waists and overlap less when worn on people with large waists. The upper edge of the waistband should be just below the navel. Fasten the outer elastic bands at your side.
3. Adjust suspenders for a comfortable fit. For maximum effectiveness, the ProFlex[®] back support must be properly positioned.
4. Before lifting, grasp the outer elastic bands and pull forward as far as possible. For model 1400UN, overlap of the right elastic band over the left elastic band is acceptable for small waist sizes. Secure bands in a comfortable position.

NOTE:
The back support should be worn loosely when not lifting.

CARE INSTRUCTIONS:
Machine wash cold; air dry.

* We cannot guarantee wearing personal protective equipment will prevent injury.

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